

Using the Ubiquity of the Cell Phone to Record Physiological Activities

Mike Sinclair¹, Aiden Doherty^{1,2}, Kristin M. Tolle¹

¹ Microsoft Research

² Dublin City University

- Background
 - Cell Phone Ubiquity
 - Lifelogging
- Cell Phone Data Logger
 - A Cell Phone Data Logging Framework
- My Physiological Diary
 - Reviewing Physiological Data Using Contextual Information
- Ongoing Work

Cell Phone Ubiquity

- 2.5+ billion Cell Phones in World (Approx 1bn PCs);
 90% of the world's population is in cell tower range
- Almost 70% of new cell phone subscriptions come from developing nations (Source: International Telecommunications Union)
- Bluetooth is now standard on most cell phones


Cell Phone as a Platform for Healthcare RFP


- 14 Projects awarded funding –
 Final reports for projects this week
- Several in posters, talks and demos at the Faculty Summit
- Culminating in an mHealth event taking place in Washington, DC in October



How Often Do You Visit Your Doctor?
























Project Aims

- Utilize cell phone ubiquity
 - Logging platform on Windows Mobile devices
 - Framework allows easy integration of new BT sensors
- Reviewing physiological values
 - Interface to monitor, analyze & browse through huge volumes of sensor data
 - "Individualize" medical baselines

SmartLogger Overview







Images



Cellular Wi-Fi LAN ActiveSync







Body Temperature

Easily Include New Sensors

Research





Location



Images







*

*

Sleep apnea neck cuff with oximeter, accelerometer, microphone, & gsr sensors



Zephyr HxM: An example of a modern biometric sensor sending out BT readings

How to Review Lots of Data?

- Physiological data:
- Little emphasis on visualization




IEEE TRANSACTIONS ON INFORMATION TECHNOLOGY IN BIOMEDICINE, VOL. 8, NO. 4, DECEMBER 2004 439 A Wireless PDA-Based Physiological Monitoring System for Patient Transport Yuan-Hsiang Lin, I-Chien Jan, Patrick Chow-In Ko, Yen-Yu Chen, Jau-Min Wong, and Gwo-Jen Jan


Our Take...

- To effectively help people understand their physiological data it must be:
 - Passively logged cell phone (no manual data entry)
 - Data gives potential associative "cues" (context+images)
 - Queryable on "temporal" axes (calendar constraints)
 - Highlight more "distinctive" events (charts)


My Physiological Diary




My Physiological Diary: Location Context




My Physiological Diary: Image Context




My Physiological Diary: Compare Across Days/Months/Years




My Physiological Diary: Display Normalised Values



My Physiological Diary: Delve Deeper into Data



My Physiological Diary: Adaptively Query Based on Time



Results

- Allows sensor device researchers to concentrate on their hardware/chemistry/physics strengths
- Will allow machine learning researchers to easily aggregate data to apply their techniques
- Will allow health conscious individuals to more easily make sense of the data they've been collecting

Use Case – Sleep Apnea

- 12 million people in USA have sleep apnea
- Process of diagnosis can involve going to "sleep lab"
- In preliminary discussions with Sleep Disorders Center in UW Medical School

Device Details

- Sleep apnea neck cuff
- Sensors:
 - Pulse oximeter, 3D accelerometer, microphone (breath sounds and pulse), & galvanized skin response
 - Considering EEG
 - All Bluetooth enabled






Future challenge – Security

Body Temperature

Images



Future – Health Vault





Location



Images













Future – Symptom detection



Future – Zigbee











Images







Data Store





Conclusions

- Utilising cell phone ubiquity
 - Logging platform on Windows Mobile devices
 - Framework allows easy integration of new BT sensors
- Reviewing physiological values
 - Interface to monitor, analyse & browse through huge volumes of sensor data
 - "Individualize" medical baselines

Lots of exciting future directions!!!



Research

Microsoft®

Your potential. Our passion.™

© 2009 Microsoft Corporation. All rights reserved. Microsoft, Windows, Windows Vista and other product names are or may be registered trademarks and/or trademarks in the U.S. and/or other countries. The information herein is for informational purposes only and represents the current view of Microsoft Corporation as of the date of this presentation. Because Microsoft must respond to changing market conditions, it should not be interpreted to be a commitment on the part of Microsoft, and Microsoft cannot guarantee the accuracy of any information provided after the date of this presentation. MICROSOFT MAKES NO WARRANTIES, EXPRESS, IMPLIED OR STATUTORY, AS TO THE INFORMATION IN THIS PRESENTATION.